

# INSTRUCTIONS FOR USE ANALGESIC DEMAND VALVE WITH INTEGRAL REGULATOR















	<b>PAGE</b>
1. SYMBOLS	2
2. INTENDED USE OF DEVICE	2
3. TECHNICAL DESCRIPTION	2
4. WARNINGS	4
5. INSTRUCTIONS FOR USE	5&6
6. CLEANING AND DISINFECTION	7
7. ACCESSORIES	7
8. INSPECTION AND USER CHECKS	7
9. MAINTENANCE	8
10. SERIAL NUMBER	8
11. INTENDED LIFE	9
12. DEMAND VALVE SPECIFICATION	9
13. REGULATOR PERFORMANCE	9
14. SPARE PARTS	9
15. APPLICABLE STANDARDS	10
16. WARRANTY	10
17. MANUFACTURER AND AUTHORISED REPRESENTATIVE DETAILS	12

**READ THESE INSTRUCTIONS BEFORE USING THE EQUIPMENT**



## 1. SYMBOLS

	Indicates a potentially hazardous situation which could result in injury to the user or others, if not avoided.
	Read these instructions before using the equipment
	Do not use any form of grease or oil
	No Smoking
	Do not use this device near any source of ignition
	Product part number
	Serial Number
	Medical device
	Date of manufacture
	Authorised Representative details
	Manufactured by
	Next Maintenance or Service Due date


## 2. INTENDED USE OF DEVICE

The B.N.O.S. Meditech Analgesic Demand Valve with Integral Regulator is intended for the self-administration of N<sub>2</sub>O/O<sub>2</sub> 50/50% V/V, an analgesic gas, sold in the UK under the trade mark of Entonox®. N<sub>2</sub>O/O<sub>2</sub> should only be used for short-term/temporary pain relief by self-administration of the demand valve user and under supervision of suitably qualified/trained medical personnel. Please refer to any rules or restrictions that may cover the use of N<sub>2</sub>O/O<sub>2</sub> in the specific country of use of the device.

The B.N.O.S. Meditech Demand Valve with Integral Regulator is designed for use with either reusable (autoclavable) patient inhalation masks or single use patient inhalation masks/mouthpieces. It is highly recommended that a single use filter is used between the demand valve and the facemask/mouthpiece to minimise the risk of contamination of the device and subsequent cross contamination between users. Please ensure that single use filters/facemasks/mouthpieces are correctly specified for use with N<sub>2</sub>O/O<sub>2</sub> 50/50% V/V. Filters/facemasks/mouthpieces require a 22mm female connection complying with BS EN ISO 5356-1:2015.

## 3. TECHNICAL DESCRIPTION

The principle of operation of the demand valve is similar to the original design of equipment used in the diving market. A diaphragm (a) contained within the demand valve housing (b) acts upon a tilt valve assembly (c). When the user exerts a negative pressure at the patient connection port (d) by inhaling, the diaphragm actuates the tilt valve due to the negative pressure created within the



housing subsequently supplying a flow of gas at the patient connection of the device. (See Fig 1) The greater the negative pressure exerted on the device, the higher the flow rate available at the patient connection port. The flow rate is entirely related to the inspiratory effort of the user.

The demand valve has been specifically designed to be used with low cost single use generic mouth piece and filter systems that are widely available using standard connections that conform to BS EN ISO 5356-1.

An Accessory (1), the AGSS Fitting Part No.673-9023-00 is available to attach to the demand valve for use with Anaesthetic Gas Scavenging Systems. (See Fig 2)

This Integrated N<sub>2</sub>O/O<sub>2</sub> Regulator is designed for use with N<sub>2</sub>O/O<sub>2</sub> cylinders fitted with a pin index valve of the correct offset pattern for medical N<sub>2</sub>O/O<sub>2</sub>.


Figure A: Integral Pin Indexed Regulator.

- a. Pressure Gauge
- b. Regulator Body
- c. Yoke sealing washer (Bodok seal)
- d. N<sub>2</sub>O/O<sub>2</sub> Gas Specific Pin Arrangement
- e. Yoke
- f. Yoke Screw
- g. Proprietary Fitting
- h. N<sub>2</sub>O/O<sub>2</sub> Hose






Figure B: Integral Pin Indexed Regulator assembled to N<sub>2</sub>O/O<sub>2</sub> Cylinder

- a. Pressure Gauge
- b. Yoke
- c. Yoke Screw
- d. Medical Gas cylinder valve body
- e. Medical Gas cylinder
- f. Medical Gas cylinder valve
- g. Proprietary Fitting
- h. N<sub>2</sub>O/O<sub>2</sub> Hose



## 4. WARNINGS!

- 4.1 Please read these instructions before use. Do not use the demand valve if you do not understand the instructions detailed in this user guide.
- 4.2  Do not use any form of grease or oil (hydrocarbon based substances) with this device.
- 4.3  Do not use this device near any source of ignition e.g., naked flame, electrically powered heaters, cigarettes etc.
- 4.4  Do not smoke any products including tobacco when using this device.
- 4.5 N<sub>2</sub>O/O<sub>2</sub> 50/50% V/V is a medicinal gas and its use must be supervised by correctly trained medical staff and/or prescribed for use where applicable by a suitably authorised medically trained individual e.g., Doctor.
- 4.6 Do not use the “Press to Flush/Test” button when the device is in use.
- 4.7 Continuous exposure to Nitrous Oxide can be harmful. Where applicable consult available standards relating to exposure levels to Nitrous Oxide. A gas scavenging accessory is available for the demand valve (AGSS fitting) for use in environments where adequate ventilation is not available i.e., indoors, inside an ambulance. The AGSS should be used where applicable.
- 4.8 Ensure sufficient N<sub>2</sub>O/O<sub>2</sub> 50/50% V/V is available for use before using the product.
- 4.9 If using medical gas cylinders, ensure that they are adequately secured or stowed at all times.
- 4.10 Once Demand Valve with integral regulator is correctly connected to the cylinder, always ensure that the cylinder valve is opened slowly.
- 4.11 The device must not be used with harness or retained face mask.
- 4.12 When the demand valve is not in use, turn off the regulator where appropriate and disconnect the valve from the medical gas supply.
- 4.13 Check the demand valve and supply hose regularly for leaks. Remove any leaking device from service immediately.
- 4.14 Ensure that the hose is always arranged in such a way that it cannot be damaged and does not cause any form of obstruction or hazard.
- 4.15 Do not use the device in conjunction with a gas supply from a cylinder which is at sub-zero temperatures, because the gas mixture in the cylinder will separate at these low temperatures.
- 4.16 Important information regarding the storage of cylinders containing pre-mixed cylinders of oxygen and nitrous oxide.
- 4.17 For Domiciliary use (up to and including 5.5 L water capacity). Before use, ensure that the cylinder has been adequately warmed by retaining it in a warm place (above 10°C) for at least 2h, or by placing it in warm water at body temperature for 5min. Then agitate the contents by inverting the cylinder three times. Never use hot water. Avoid wetting the valve.
- 4.18 For Hospital use. On receipt, the cylinders shall be date marked and stored horizontally for 24h in a room maintained at a temperature above 10°C, but not exceeding 38°C, before being used. During delivery from the store room to the final destination in the hospital, it is essential that the cylinder shall not be exposed to freezing temperatures for more than 2min to 3min. If exposure to freezing temperatures occurs for any longer period, carry out the procedure in k)<sup>1)</sup> of BS 4272-2:1996. Special instructions for hospitals shall be displayed where O<sub>2</sub> and N<sub>2</sub>O mixtures are stored and the attention of all personnel handling the cylinders should be directed to them

## 5. INSTRUCTIONS FOR USE

### 5.1 ATTACHING TO A CYLINDER

- 5.1.1 Inspect the Pin Indexed Regulator (Refer to Figure: A) and check that the Yoke sealing washer (c) is in place and in serviceable condition.
- 5.1.2 Check your N<sub>2</sub>O/O<sub>2</sub> Cylinder to ensure that it is correctly labelled as medical nitrous oxide and oxygen. Ensure that any new gas cylinders are fitted with a tamper evident seal. Remove all traces of the seal. The N<sub>2</sub>O/O<sub>2</sub> Gas Specific Pin Arrangement (referring to Figure: A (d) prevents incorrect attachment of the wrong type of gas cylinder and will locate only into the matching hole on the face of the N<sub>2</sub>O/O<sub>2</sub> cylinder valve.
- 5.1.3 Referring to Figure B: Position the Regulator Yoke (b) over the Cylinder valve body (d), having first loosened the Yoke Screw (c). Position the N<sub>2</sub>O/O<sub>2</sub> Gas Specific Pin Arrangement (a) into the locating hole and tighten Yoke Screw (c) securely, ensuring that the Pressure Gauge (a) is in the correct orientation. We recommend the orientation shown in Figure B where the hose (h) is hanging downwards. During assembly ensure the regulator gauge is not compressed against the cylinder valve, which can damage the gauge and prevent the Yoke sealing washer functioning. When the assembled regulator and cylinder is placed in the kit bag, the cylinder gauge can be easily seen and the cylinder valve (f) can be operated.
- 5.1.4 Referring to Figure B: Open the Cylinder Valve (f) slowly (i.e. facing away from the operator or other personnel). If gas escapes around the Yoke sealing washer, turn off the cylinder at the Cylinder Valve (f) and tighten the Yoke Screw (c). Again, slowly open the cylinder valve, if no gas leak is heard, open the Cylinder Valve one full turn and check that the Pressure Gauge (a) shows that the cylinder is full.
- 5.1.5 If a leak of gas is still heard, turn off the Cylinder Valve and remove the Regulator. Remove and replace the Yoke sealing washer (Refer to Figure: A) (c) and repeat from step 5.1.1.
- 5.1.6 After use the cylinder valve should be turned off.


### 5.2 TO CHANGE A USED CYLINDER

- 5.2.1 Turn off the Cylinder Valve (Fig: B) by the Cylinder Valve Control (f)
- 5.2.2 De-pressurise the regulator by pressing the "Flush/Test Button on the Demand Valve". The pressure gauge (a) reading will move to zero. Unscrew the Regulator Yoke Screw (c) and remove the empty Cylinder ensuring the Yoke sealing washer (Fig: A) (c) remains in place on a regulator.
- 5.2.3 Repeat the instructions for "Attaching to a cylinder". (Section 5.1)

### 5.3 TO USE DEMAND VALVE

- 5.3.1 Before connecting to a pressurised gas source, please ensure unit is clean and in good condition. If you have any doubts about the condition, please do not connect or use the device. Please check that the inhalation diaphragm is fitted to the device and seated correctly (lying flat against the body of the device).
- 5.3.2 Ensure that the supply hose attached to the demand valve is connected to an N<sub>2</sub>O/O<sub>2</sub> Medical Gas supply that uses the correct respective Medical Gas specific fitting.
- 5.3.3 If using a N<sub>2</sub>O/O<sub>2</sub> regulator, please check that the cylinder has adequate gas and that the cylinder valve is turned on. It is suggested that the portable cylinder is at least half full and that additional cylinders are available if required.

- 5.3.4 When using an AGSS, ensure that the AGSS fitting is correctly attached and seated on the demand valve and that the scavenging pipe is correctly and firmly attached.
- 5.3.5 Before administering N<sub>2</sub>O/O<sub>2</sub> to the patient, please check the availability of gas by briefly pressing the 'Flush/Test' button, located on the top of the demand valve. Please check that there are no audible sounds of gas flow once the 'Flush/Test' button has been released.
- 5.3.6 It is strongly advised to use a single use filter with the demand valve to minimise the risk of contamination of the device and subsequent cross contamination between users. Connect a universal patient inhalation mask or single use patient inhalation mask/mouthpiece to the single use filter.
- 5.3.7 The demand valve is now ready for use. The user should be instructed to inhale through the mouthpiece or face mask connected to the patient connection port fitted to the demand valve. Doing so will create a negative pressure on the patient connection port which will activate the supply of Analgesic Gas. The larger the breath taken by the user, the greater the supply of N<sub>2</sub>O/O<sub>2</sub>. Please advise the user not to press the 'Flush/Test' button whilst the device is in use.
- 5.3.8 The user is able to exhale through the mouthpiece or face mask.
- 5.3.9 The demand valve will stop supplying gas when the user stops inhaling. The exhaled gas will exit through the demand valve via the exhalation port.
- 5.3.10 The device should continue to be used under the supervision of medically trained personnel.
- 5.3.11 Stop use of the demand valve if the suitably trained supervisor or user detects a change in performance of the device.
- 5.3.12 The duration of use of the device is to be determined by the medically trained professionals who are supervising the use of the demand valve.



## 6. CLEANING AND DISINFECTION

- 6.1 Ensure demand valve with integral regulator is disconnected from the gas supply before cleaning.
- 6.2 Clean the device before first use and then subsequently after every use, ensure the demand valve, including the handle is wiped over thoroughly with a disinfecting wipe. For disinfection purposes a chlorine dioxide-based product (e.g., the Tristel Wipes and solution system) can be used, at a nominal concentration of 0.02% wt/vol. The concentration refers to chlorine dioxide in water. Follow the manufacturer's directions for safe use.
- 6.3 The demand valve with integral regulator is not suitable for autoclaving.
- 6.4 The demand valve with integral regulator should be thoroughly dried before storage.

## 7. ACCESSORIES

Part No.	673-9023-00	Description: AGSS Fitting
----------	-------------	---------------------------

## 8. INSPECTION AND USER CHECKS

### 8.1. INSPECTION

- 8.1.1 The regulator should be inspected for damaged or broken components and contamination after each use. Check that the Pressure Gauge (a) and Proprietary Fitting (g) are all tight and undamaged. Check that any controls operate correctly and that the Yoke Screw Fig A.(f) operates smoothly.
- 8.1.2 If contaminated, the Regulator should be cleaned in accordance with Section 6.
- 8.1.3 If damaged, the Regulator should be withdrawn from service and returned to B.N.O.S. Meditech for assessment and repair.

### 8.2. USER CHECKS

- 8.2.1 Connect the Demand Valve with Integral Regulator to a suitable Medical Gas cylinder. Check that the Pressure Gauge (a) reads zero. Slowly open the cylinder valve. Check that the Pressure Gauge indicates that there is pressure in the cylinder. If it continues to read zero, use a cylinder known to contain Medical Gas at a minimum pressure of 50 bar,
- 8.2.2 Check for audible leaks.
- 8.2.3 Close the cylinder valve.
- 8.2.4 If the Regulator fails the check, it should be withdrawn from service and returned to B.N.O.S. Meditech Ltd. for assessment and repair.

## 9. MAINTENANCE

- 9.1 Maintenance must be carried out on the unit on a five-yearly basis by B.N.O.S. Meditech or engineers certified by B.N.O.S. Meditech. This activity involves dismantling the unit and replacing all internal seals, diaphragm, medical gas supply hose and any components which show significant wear and tear.
- 9.2 Performance should also be checked on a five-yearly basis, using suitable test equipment. The negative pressure required to generate specific flows from the device must be tested. If required, B.N.O.S. Meditech Ltd. can advise on suitable test equipment.
- 9.3 For demand valves in frequent use, for example, an Ambulance Service or a Maternity Unit, personnel trained to supervise the use of a demand valve should refer to 5.3.5 of this Instruction For Use document and carry out the check as described on a regular basis. For less frequent users, i.e. industrial safety use, a leak test should be periodically performed by applying an oxygen compatible leak test solution to all medical gas supply hose fittings and joints, where applicable. The suggested period is annually.
- 9.4 The medical gas supply hose MUST be changed every five years. Please contact the manufacturer to determine the age of the hose fitted to the demand valve.
- 9.5 B.N.O.S. Meditech Ltd. offers training and certification on the service, repair and preventative maintenance of its products.

## 10. SERIAL NUMBER

The Meditech Demand Valve with Integral Regulator has 2 serial numbers and 1 overall product part number. The Demand Valve has its own serial number, as does the integral regulator.

The Demand Valve Serial No. and the overall product code is to be found on the handle label of the demand valve. It consists of four sections, the first letters expressing the general type of device, followed by numbers for the month and year of manufacture and lastly a set of up to five numbers representing the individual “number” of the unit and differentiating units built in the same month.

Demand Valve Example: SN EDV121700001

EDV	12	17	00001
Entonox® Demand Valve	Month of Manufacture	Year of Manufacture	Number of Unit

The serial number for the integral regulator is to be found on the “barrel” label of the regulator.

Integral Regulator Example: SN R121800001

R	03	18	00001
Regulator	Month of Manufacture	Year of Manufacture	Number of Unit



## 11. INTENDED LIFE

The demand valve with integral regulator is designed to have a product life span of 10 years, excluding abuse and/or damage to the device.

## 12. DEMAND VALVE SPECIFICATION

Inspiratory Resistance <1.5 kPa @ 200 l/min  
<0.25 kPa @ 10 l/min

Expiratory Resistance <0.1kPa @ 12 l/min  
<0.6kPa @ 120 l/min

Transport and Storage Temperature: -20°C to 60°C

Operating Temperature: 5°C to 40°C

Supply Pressure @ 100l/min Nominal 400 kPa  
Maximum 600 kPa  
Minimum 280 kPa

## 13. REGULATOR PERFORMANCE

Rated maximum inlet pressure (P1): 200 bar

Rated outlet pressure (P2): 4 bar

Test inlet pressure (P3): 8 bar

Standard discharge (Q1): 50 lpm (litres per minute)

Variation of outlet pressure (P2) when inlet pressure is varied from P1 to P3 at a flow of Q1: 4%

Temperature ranges: Storage: -40°C to +60°C Operating: -20°C to +60°C

## 14. SPARE PARTS

Description	Part Number
Inhalation Diaphragm	033-1030-00
Protective boot for pressure gauge	633-0027-00
Yoke sealing washer (also known as bodok seal)	033-1021-00

## 15. APPLICABLE STANDARDS

B.N.O.S. Meditech Ltd. is an ISO 13485 certified company.

### EC – DECLARATION OF CONFORMITY CE 2797

These demand valves, class IIa medical devices are supplied in accordance with the Medical Devices Directives 93/42/EEC. BNOS Meditech Oxygen Regulators are supplied in conformity under a quality system to meet Annex II of the Medical Devices Directive 93/42/EEC. Oxygen Regulators are classified as Class IIb Medical Devices.

The above quality system has been inspected by the Notified Body Ref: CE 2797 being BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

### Other Standards

The following National & International Standards apply to the device:

Standard Number:	Title:
BS 4272-2:1996	Anaesthetic and Analgesic machines. Specification for intermittent (demand) flow analgesic machines for use with 50/50% V/V Nitrous Oxide and Oxygen.
BS EN ISO 5356-1:2015	Anaesthetic and respiratory equipment, cones and connectors, cones and sockets.
BS EN ISO 15001:2011	Anaesthetic and respiratory equipment. Compatibility with Oxygen.
BS EN ISO 5359:2014+A1:2017	Low pressure hose assemblies for use with medical gases.
BS 5682:2015	Probes (quick connectors) for use with medical gas pipeline systems.
BS EN ISO 14971:2019	Medical Devices. Application of risk management to medical devices.
BS EN ISO 10524-1:2019	Pressure regulators for use with medical gases – Part 1: Pressure regulators and pressure regulators with flow-metering devices
BS EN ISO 15223-1:2021	Medical Devices. Symbols to be used with medical device labels, labelling and information to be supplied. General requirements.

It is also certified that the equipment listed above fully complies with all the required mandatory standards and the performance, specifications, standards and sources agreed and contracted for this order

#### IMPORTANT NOTICE

Manufacturer's Warranty is for a period of 5 years and includes parts and labour. It does not include transport costs. The responsibility and cost of returning and collecting the unit from the manufacturer or their authorised representative is the owners.

B.N.O.S Meditech Ltd reserves the right to change design without prior notice.

**Any disassembly of the device beyond that detailed in this manual will invalidate the warranty and the manufacturers disclaim any liability for products that have undergone unauthorised repair.**

**BLANK  
PAGE**

Entonox® is a registered trademark of The Linde Group.

**medITECH**  
MEDICAL GAS EQUIPMENT

**CE 2797**

**COMPANY CONTACT DETAILS**

This device is designed and manufactured by:

**B.N.O.S. MEDITECH LTD**

9 FIFTH AVENUE, BLUEBRIDGE IND EST  
HALSTEAD, ESSEX, CO9 2SZ, ENGLAND

Tel: +44 (0)1787 479 475

Fax: +44 (0)1787 477 747

E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)

Website: [www.meditech.uk.com](http://www.meditech.uk.com)



**EU REPRESENTATIVE:**

MEDICAL DEVICE MANAGEMENT LTD  
BLOCK B, THE CRESCENT BUILDING  
NORTHWOOD, SANTRY  
DUBLIN 9, D09 C6X8  
IRELAND