


2L Medical Oxygen Cylinder

BNOS Meditech has a Wholesalers Dealers License, allowing the sale and distribution of Medical Oxygen. Authorisation Holder Number: WDA (H) 44056



SPECIFICATION

| | |
|---------------------------------|-------------------------|
| Filling Pressure:..... | 200 BAR |
| Test Pressure:..... | 300 BAR |
| Material:..... | Aluminum Alloy 6061 |
| Diameter:..... | ~114mm |
| Weight(filled with valve):..... | ~4Kg |
| Weight(Empty with valve):..... | ~2.8Kg |
| Height:..... | ~ 444.5mm |
| Colour Code:..... | White* - EU Standard |
| Standard:..... | BS EN ISO 7866:2012 |
| Valve: | Pin Index or Bullnose** |
| Volume (L):..... | 400 Litres |
| Thread: | 17E Tapered |

* Other colours available if required

** Other valve types available if required

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Entonox is a registered Trademark of the Linde Group.

Distributor: