

# medITECH

MEDICAL GAS EQUIPMENT

## Meditech Medical Gas Hose Assemblies & Hose Adaptors



- CE & UKCA Marked
- Gas specific NIST connectors
- Phthalate free
- Non - magnetic
- Anti-microbial
- Different hose lengths available
- Anti-static hose
- 5 year service life

### AVAILABLE GAS TYPES:

• CO<sub>2</sub> • O<sub>2</sub> • MA4 • SA7 • O<sub>2</sub>/N<sub>2</sub>O • N<sub>2</sub>O • VACUUM

Create a Hose Assembly to suit your specific requirements by following the steps below to generate your **HOSE ASSEMBLY PART NUMBER**.


**STEP 1** Choose the **GAS TYPE** that you require

**STEP 2** Choose the **INLET / PROBE TYPE** you require

**STEP 3** Choose the **HOSE LENGTH** you require

**STEP 4** Choose the **OUTLET TYPE** you require

## HOSE ASSEMBLY PART NUMBERING SYSTEM



**GAS TYPES**

**O<sub>2</sub>   O<sub>2</sub>/N<sub>2</sub>O   N<sub>2</sub>O   SA7   MA4   VACUUM   CO<sub>2</sub>**

**PROBE TYPES**




DISS NUT (F)



**HOSE LENGTH**

**1M   2M   3M   4M   5M   6M   7M**

**OUTLET TYPES**



# medITECH

## MEDICAL GAS EQUIPMENT

**TO PURCHASE OR FOR FURTHER DETAILS:**

Contact us by email or telephone for more details,  
including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**

**9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex, CO9 2SZ, England, UK.**

**TEL: +44 (0) 1787 479 475**

**FAX: +44 (0) 1787 477747**

**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**

**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**

Meditech Hose Assemblies are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.



Distributor: