

# mediTECH

MEDICAL GAS EQUIPMENT

## Meditech Medical Gas Hose Assemblies & Hose Adaptors



- CE & UKCA Marked
- Phthalate free
- Anti-microbial
- Anti-static hose
- Gas specific NIST connectors
- Non - magnetic
- Different hose lengths available
- 5 year service life

### AVAILABLE GAS TYPES:

• CO<sub>2</sub> • O<sub>2</sub> • MA4 • SA7 • O<sub>2</sub>/N<sub>2</sub>O • N<sub>2</sub>O • VACUUM

Create a Hose Assembly to suit your specific requirements by following the steps below to generate your **HOSE ASSEMBLY PART NUMBER**.

**STEP 1** Choose the **GAS TYPE** that you require

**STEP 2** Choose the **INLET / PROBE TYPE** you require

**STEP 3** Choose the **HOSE LENGTH** you require

**STEP 4** Choose the **OUTLET TYPE** you require

## HOSE ASSEMBLY PART NUMBERING SYSTEM

### Gas type

151 = O<sub>2</sub>  
152 = O<sub>2</sub>/N<sub>2</sub>O  
153 = N<sub>2</sub>O  
154 = SA7  
155 = MA4  
156 = Vacuum  
157 = CO<sub>2</sub>

### Hose Length

50 = 0.5m      04 = 5m  
75 = 0.75m    05 = 6m  
00 = 1m        06 = 7m  
51 = 1.5m      07 = 8m  
01 = 2m  
02 = 3m  
03 = 4m



### Probe Type

|                             |                        |
|-----------------------------|------------------------|
| 00 = BS Probe               | 11 = Ohmeda            |
| 01 = BS Probe (90)          | 12 = Puritan           |
| 02 = NIST Probe + Nut       | 13 = Chemetron         |
| 03 = NIST Probe + Nut (90°) | 14 = DISS NUT (F)      |
| 04 = Afnor                  | 15 = DISS (90°)NUT (F) |
| 05 = AGA                    | 16 = 3/8" BSP Nut      |
| 06 = Czech                  | 17 = DISS (M)          |
| 07 = SIS                    | 18 = DIN               |
| 08 = DIN (90)               | 19 = Mini Schrader     |
| 09 = Unifor                 | 20 = BS Probe - SRV    |
| 10 = DISS HAND (F)          | 21 = RA AGA            |






### Outlet Type

|                             |                         |
|-----------------------------|-------------------------|
| 00 = BS Schrader            | 15 = 3/8" BSP Nut (90°) |
| 01 = NIST Nut + Probe       | 16 = BS Probe           |
| 02 = NIST Probe + Nut (90°) | 17 = DISS HAND (F)      |
| 03 = Afnor                  | 18 = DISS (M)           |
| 04 = AGA                    | 19 = 3/8" BSP           |
| 05 = Czech                  | 20 = Afnor (90°)        |
| 06 = SIS                    | 21 = OHIO               |
| 07 = DIN                    | 22 = 1/8" BSP (M)       |
| 08 = Unifor                 | 23 = 1/4" BSP (F)       |
| 09 = DISS NUT (F)           | 24 = 1/4" BSP (M)       |
| 10 = Ohmeda                 | 25 = J512 FML Nut       |
| 11 = Chemetron              |                         |
| 12 = Puritan                |                         |
| 13 = DISS (90°)(MALE)       |                         |
| 14 = 3/8" BSP Nut           |                         |

GAS TYPES

O<sub>2</sub>   O<sub>2</sub>/N<sub>2</sub>O   N<sub>2</sub>O   SA7   MA4   VACUUM   CO<sub>2</sub>




PROBE TYPES

|                                                                                     |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|    |  |  |  |  |   |  |  |
| BS Probe                                                                            | BS Probe (SRV)                                                                    | BS Probe (90°)                                                                    | NIST Probe + Nut                                                                  | NIST Probe + Nut (90°)                                                            | AFNOR                                                                               | AGA                                                                                 | Czech                                                                               |
|    |  |                                                                                   |                                                                                   |  |  |  |  |
|    |  |  |  |  |   |  |  |
| SISS                                                                                | DIN                                                                               | UNIFOR                                                                            | DISS HAND (F)                                                                     | Ohmeda                                                                            | Puritan                                                                             | Chemetron                                                                           | DISS Male                                                                           |
|                                                                                     |  |  |  |  |  |  |  |
|   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
| DISS NUT (F)                                                                        |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|  |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |

HOSE LENGTH

1M   2M   3M   4M   5M   6M   7M

OUTLET TYPES

|                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                       |                                                                                       |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |  |   |  |  |
| BS Schrader                                                                         | NIST Nut + Probe                                                                    | NIST Probe + Nut (90°)                                                              | AFNOR                                                                               | DIN                                                                                 | Czech                                                                                 | AGA                                                                                   | SIS                                                                                   |
|  |                                                                                     |                                                                                     |  |  |  |  |  |
|  |  |  |  |  |   |  |  |
| UNIFOR                                                                              | DISS (M)                                                                            | Ohmeda                                                                              | Chemetron                                                                           | Puritan                                                                             | DISS HAND (F)                                                                         | DISS NUT (F)                                                                          | OHIO                                                                                  |
|  |  |  |  |  |  |  |  |



# medITECH

## MEDICAL GAS EQUIPMENT

**TO PURCHASE OR FOR FURTHER DETAILS:**

Contact us by email or telephone for more details,  
including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**

**9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex, CO9 2SZ, England, UK.**

**TEL: +44 (0) 1787 479 475**

**FAX: +44 (0) 1787 477747**

**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**

**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**



Meditech Hose Assemblies are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor: