

medITECH

MEDICAL GAS EQUIPMENT

Meditech Medical Gas Hose Assemblies & Hose Adaptors



- CE & UKCA Marked
- Phthalate Free
- Anti-microbial
- Anti-static Hose
- Gas Specific NIST connectors
- Non - magnetic
- Different hose lengths available
- 5 year service life

AVAILABLE GAS TYPES:

• O₂ • MA4 • SA7 • O₂/N₂O • N₂O • VACUUM

Create a Hose Assembly to suit your specific requirements by following the steps below to generate your **HOSE ASSEMBLY PART NUMBER**.


STEP 1 Choose the **GAS TYPE** that you require

STEP 2 Choose the **INLET / PROBE TYPE** you require

STEP 3 Choose the **HOSE LENGTH** you require

STEP 4 Choose the **OUTLET TYPE** you require































HOSE ASSEMBLY PART NUMBERING SYSTEM



GAS TYPES

O₂ O₂/N₂O N₂O MA7 MA4 VACUUM CO₂































PROBE TYPES


							
BS Probe	BS Probe (90°)	NIST Probe + Nut	NIST Probe + Nut (90°)	AFNOR	AGA	Czech	SIS
							
							
DIN	UNIFOR	DISS HAND (F)	Ohmeda	Puritan	Chemetron	DISS Male	DISS NUT (F)
							

HOSE LENGTH

1M 2M 3M 4M 5M 6M 7M

OUTLET TYPES

							
BS Schrader	NIST Nut + Probe	NIST Probe + Nut (90°)	AFNOR	DIN	Czech	AGA	SIS
							
							
UNIFOR	DISS (M)	Ohmeda	Chemetron	Puritan	DISS HAND (F)	DISS NUT (F)	OHIO
							



medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.**

TEL: +44 (0) 1787 479 475

FAX: +44 (0) 1787 477747

EMAIL: sales@meditech.uk.com

WEBSITE: www.meditech.uk.com

Meditech Hose Assemblies are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked accordance with the UK Medical Devices Regulations 2002. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Entonox is a registered Trademark of the Linde Group.



Distributor: