


2L Entonox® Cylinder

BNOS Meditech has a Wholesalers Dealers License, allowing the sale and distribution of Medical Gas. Our entonox® cylinders are available either filled or empty with a pin index valve. Other valves are available upon request.

Authorisation Holder Number: WDA (H) 44056



SPECIFICATION

Filling Pressure:.....	137 BAR
Material:.....	Aluminum Alloy
Diameter:.....	~111mm
Weight(filled with valve):.....	~4Kg
Weight(Empty with valve):.....	~2.8Kg
Height:.....	~ 420mm
Colour Code:.....	White with Blue Shoulder - EU Standard
Standard:.....	BS EN ISO 7866:2012
Valve:	Pin Index
Volume (L):.....	274 Litres
Thread:	17E Tapered

meditech

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.