

medITECH

MEDICAL GAS EQUIPMENT

50% Nitrous Oxide / 50% Medical Oxygen 2L Cylinder

BNOS Meditech has a Wholesalers Dealers License, allowing the sale and distribution of Medical Gas. Our cylinders are available with different valve types, either filled or empty.

**Wholesale Authorisation Holder Number:
WDA (H) 44056**



SPECIFICATION

Filling Pressure:.....	137 BAR
Material:.....	Aluminum Alloy
Diameter:.....	~102 mm
Weight(filled with valve):.....	~4.2Kg
Weight(Empty with valve):.....	~3Kg
Height:.....	~525 mm
Colour Code:.....	White with Blue Shoulder - EU Standard
Standard:	BS EN ISO 7866:2012
Valve:	Integral
Volume (L):.....	428 Litres
Outlet Type:.....	Therapy flow (0-15lpm), 4 bar BS Schrader

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.