



MicroVENT™ EUROPEAN RESUSCITATOR

- **Pneumatic Oxygen powered resuscitator.**
- **Automatic (time cycled) and manual trigger.**
- **Suitable for adults & children.**
- **Easy to use.**
- **Lightweight.**
- **Compact.**

A hand-held, lightweight resuscitator with compact dimensions, the Microvent European Resuscitator has ergonomic styling to give easy operation with controls in a natural position. It is constructed with a tough reinforced plastic casing to withstand use in difficult situations.

The Microvent patient valve features 22/15mm connections to take standard airway devices and face masks. The oxygen supply hose is 2 metres long and made of resilient reinforced plastics to retain flexibility in cold conditions.

A British Standard (BS) medical oxygen supply probe enables use with all BS standard portable and piped oxygen supplies. Other international hose fittings are also available.

SPECIFICATION

MICROVENT European RESUSCITATOR

| | |
|--|---------------------------|
| Patient population range: | Adult / child above 14kg |
| Automatic operation: | Time cycled. Gas powered. |
| Automatic flow rate (L/min): | 15.5 - 21.5 |
| Automatic tidal volume (L): | 0.6 to 0.2 |
| Automatic oxygen concentration V/V:..... | 100% |
| Automatic I:E ratio:..... | 1:2 |
| Automatic frequency (per minute): | 12 - 25 |
| Manual flow rate (L/MIN): | 40 |
| Pressure relief valve with audible warning limits maximum attainable delivery pressure (KrPa):..... | 4.5 (6.0 on request) |
| Expiratory resistance (KrPa):..... | <0.5 |
| Inspiratory resistance (KrPa):..... | <0.5 |
| Resuscitator weight excluding supply hose (g):..... | 214g |
| Maximum resuscitator dimensions excluding supply hose (mm):..... | 120 x 55 x 100 |

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

Meditech microVENTs are CE marked in accordance with the Medical Devices Directive in the EU and UKCA to the Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

medITECH
MEDICAL GAS EQUIPMENT

Distributor: