

## Straight Bar Flow Selector BS Schrader Outlet

Direct BS Probe - 100mm



- **No Calibration required EVER**
- **Lightweight construction**
- **Available with International Fittings**
- **7 Year Warranty**
- **7 Year service Interval**

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Twin Flow Selector is designed to be plugged into Medical Oxygen Wall Outlets.

The Meditech Straight bar Flow Selector, BS Schrader Outlet has been designed to deliver precision flow rates, without the need for any calibration. It is fitted with a BS Schrader 4 BAR Outlet.

Can be supplied in either of two formats: Push-Fit 'Firtree' Fitting for direct attachment on a Oxygen Therapy delivery system With a threaded outlet fitting for the attachment of a humidifier.

**AVAILABLE GAS TYPES:**  
OXYGEN | MEDICAL AIR | ENTONOX

### THERAPY OUTLET TYPES



FIR TREE



9/16<sup>th</sup> THREADED  
THERAPY CONNECTOR




3/8" BSP THREADED  
THERAPY CONNECTOR



1/4" BARB

### PROBE TYPE



BS PROBE



DISS



DIN



AGA



AFNOR

SIS



CZECH



UNIFOR



### OUTLET TYPES



BS SCHRADER



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNIFOR



OHIO



# SPECIFICATION

|                                    |   |
|------------------------------------|---|
| Intended Gas: .....                | Medical Oxygen (O <sub>2</sub> , MA4, O <sub>2</sub> /N <sub>2</sub> O)   |
| Inlet pressure: .....              | 4 bar (400 kPa)   |
| Set Flow Rates in lpm: .....       | Three flow rates options available:<br>1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm<br>1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm<br>10, 20, 30, 35, 40, 50, 55, 60, 70 lpm |
| Inlet Connections: .....           | Schrader (BS 5682:2015) and other international connections   |
| Outlet (Patient) Connections. .... | Fintree style barbed outlet, 9/16 threaded outlet, 3/8" BSP threaded outlet & 1/4" barb   |
| Intended Life: .....               | 10 Years  |
| Service Interval: .....            | 7 Years   |
| Warranty Period: .....             | 7 Years   |
| Complies With: .....               | ISO 15002: 2023   |

**TO PURCHASE OR FOR FURTHER DETAILS:**  
**Contact us by email or telephone for more details,**  
**including specific part numbers, prices, etc.**

B.N.O.S. Meditech Ltd.  
9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex CO9 2SZ, England, UK.  
Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747  
E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com) • [www.meditech.uk.com](http://www.meditech.uk.com)



**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor:

Meditech Flow Selectors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Doc 784M Issue D 28/08/2024