

## Straight Bar, Rail Mounted Flow Selector 2M Supply Hose, BS Probe -100mm



- **No Calibration required EVER**
- **Lightweight construction**
- **Available with International Fittings**
- **7 Year Warranty**
- **7 Year service Interval**

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Straight Bar Flow Selector and Outlet is designed to be plugged into Medical Oxygen Wall Outlets, whilst being mounted on a medical rail.

The Meditech Straight Bar Flow Selector and Outlet has been designed to deliver precision flow rates, without the need for any calibration. It's fitted with a universal size rail mount for most rail systems. It has a standard 2M hose length (other hose lengths available upon request).

Can be supplied in either of two formats: Push-Fit 'Firtree' Fitting for direct attachment on a Oxygen Therapy delivery system With a threaded outlet fitting for the attachment of a humidifier.

### THERAPY OUTLET TYPES



FIR TREE



9/16<sup>th</sup> THREADED  
THERAPY CONNECTOR



3/8" BSP THREADED  
THERAPY CONNECTOR



1/4" BARB


### AVAILABLE GAS TYPES:

OXYGEN | MEDICAL AIR | ENTONOX

### PROBE TYPE



BS PROBE



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNIFOR



### OUTLET TYPES



BS SCHRADER



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNIFOR



OHIO



# SPECIFICATION

Intended Gas: .....	Medical Oxygen (O <sub>2</sub> , MA4, O <sub>2</sub> /N <sub>2</sub> O)
Inlet pressure: .....	4 bar (400 kPa)
Set Flow Rates in lpm: .....	Three flow rates options available: 1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm 10, 20, 30, 35, 40, 50, 55, 60, 70 lpm
Inlet Connections: .....	Schrader (BS 5682:2015) and other international connections
Outlet (Patient) Connections. ....	Firtree style barbed outlet, 9/16 threaded outlet, 3/8" BSP threaded outlet & 1/4" barb
Intended Life: .....	10 Years
Service Interval: .....	7 Years
Warranty Period: .....	7 Years
Complies With: .....	ISO 15002: 2008

## TO PURCHASE OR FOR FURTHER DETAILS:

**Contact us by email or telephone for more details, including specific part numbers, prices, etc.**

B.N.O.S. Meditech Ltd.  
9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex CO9 2SZ, England, UK.  
Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747  
E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com) • [www.meditech.uk.com](http://www.meditech.uk.com)

Meditech Flow Selectors are CE marked in accordance with the Medical Devices Directive in the EU and UKCA to the Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website. Entonox® is a trade mark of the Linde Group

**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor: