



0-3LPM

Low Flow Selector

- **No Calibration required EVER**
- **7 Year Warranty**
- **7 Year Service Interval**
- **Fir Tree Styled Barbed Outlet or Threaded Fitting**
- **Highly Accurate Pre-Set Flow Rates**

Excellent performance and robustness, coupled with very accurate pre-set flow rates.

For use in all Hospital and Emergency environments, the Low Flow Flow Selector is designed to be plugged into Medical Oxygen Wall Outlets or portable Oxygen cylinders.

The mediTECH Low Flow Flow Selector has been designed to deliver precision flow rates, without the need for any calibration.

Can be supplied in either two formats;

- Push-fit 'fir-tree' fitting for direct attachment on a Oxygen Therapy delivery system.
- Threaded Outlet fitting for the attachment of external devices.

0-1LPM



SPECIFICATION

Intended Gas:

Medical Oxygen (O₂)

Inlet pressure:

4 bar (400 kPa)

Set Flow Rates in lpm:

0.01, 0.02, 0.03, 0.04 0.05, 0.06, 0.07, 0.08, 0.09, 0.10 & 1lpm
0.02, 0.03, 0.05, 0.08, 0.12, 0.20, 0.30, 0.50, 0.75, 1 & 3lpm

Inlet Connections:

Schrader (BS 5682:2015) and
other international connections

Outlet (Patient) Connections:

Firtree style barbed outlet or Threaded fitting
(for use with humidifiers)

Intended Life:

10 Years

Service Interval:

7 years

Warranty Period:

7 years

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com



Distributor:

Meditech Flow Selectors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Doc 845M Issue C 29/09/2025