



## Flow Selector

With extended BS Probe

- **No Calibration required EVER**
- **Lightweight construction**
- **Available with International Fittings**
- **5 Year Warranty**
- **5 Year service Interval**

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Flow Selector is designed to be plugged into Medical Oxygen Wall Outlets. This Flow Selector has an extended BS Probe.

The Meditech Flow Selector has been designed to deliver precision flow rates, without the need for any calibration. It has a standard 2M hose length (other hose lengths available upon request).

Can be supplied in either of two formats: Push-Fit 'Firtree' Fitting for direct attachment on a Oxygen Therapy delivery system With a threaded outlet fitting for the attachment of a humidifier.

### OTHER THERAPY TYPES:



FIR TREE



9/16" THREADED THERAPY CONNECTOR



3/8" BSP THREADED THERAPY CONNECTOR



1/4" BARB

### AVAILABLE GAS TYPES:

- O<sub>2</sub>
- MA4
- Entonox

# SPECIFICATION

## Twin Flow Selector, Straight Bar

Intended Gas: .....	Medical Oxygen (O <sub>2</sub> , MA4, O <sub>2</sub> /N <sub>2</sub> O)
Inlet pressure: .....	4 bar (400 kPa)
Set Flow Rates in lpm: .....	Three flow rates options available: 1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm 10, 20, 30, 35, 40, 50, 55, 60, 70 lpm
Inlet Connections: .....	Schrader (BS 5682:2015) and other international connections
Outlet (Patient) Connections. ....	Firtree style barbed outlet, 9/16 threaded outlet, 3/8" BSP threaded outlet & 1/4" barb
Intended Life: .....	10 Years
Service Interval: .....	5 Years
Warranty Period: .....	5 Years
Complies With: .....	ISO 15002: 2008
Extended Probe Length: .....	7CM

### TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

#### **B.N.O.S. Meditech Ltd.**

**9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex, CO9 2SZ, England, UK.**

**TEL: +44 (0) 1787 479 475**

**FAX: +44 (0) 1787 477747**

**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**

**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**

Meditech Flow Selectors are Class IIa medical devices, CE marked in accordance with the Medical Devices Directive in the EU. All items in this range are sold with a return to manufacturer warranty against manufacturing defects.

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website. Entonox is a registered Trademark of the Linde Group.

**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor:

**CE 2797**

Doc 769M Issue A 21/04/21