

Twin Flow Selector Wall Mounted

With 2M supply hose ***Also available with Rail Mount***

- No Calibration required EVER
- Lightweight construction
- Available with International Fittings
- 5 Year Warranty • 5 Year service Interval




Excellent performance and robustness, coupled with a very lightweight design. For use in all Hospital and Emergency environments, the Twin Flow Selector is designed to be plugged into Medical Oxygen Wall Outlets, whilst being attached directly to the wall, or attached to a rail. The Meditech Twin Flow Selector has been designed to deliver precision flow rates, without the need for any calibration. Its fitted with either a universal wall mounting plate or a universal size rail mount for most rail systems. It has a standard 2M hose length. *(other hose lengths available upon request).*


Can be supplied in either of two formats:

- Push-Fit 'Firtree' Fitting for direct attachment on a Oxygen Therapy delivery system
- With a threaded outlet fitting for the attachment of a humidifier. *(Humidifiers to be purchased seperately)*

PROBE TYPES:



THERAPY TYPES:



AVAILABLE GAS TYPES:

- O2
- MA4
- Entonox

SPECIFICATION

Intended Gas:	Medical Oxygen (O ₂ , MA4, O ₂ /N ₂ O)
Inlet pressure:	4 bar (400 kPa)
Set Flow Rates in lpm:	Three flow rates options available: 1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm 10, 20, 30, 35, 40, 50, 55, 60, 70 lpm
Inlet Connections:	Schrader (BS 5682:2015) and other international connections
Outlet (Patient) Connections.	Firtree style barbed outlet, 9/16 threaded outlet, 3/8" BSP threaded outlet & 1/4" barb
Intended Life:	10 Years
Service Interval:	5 Years
Warranty Period:	5 Years
Complies With:	ISO 15002: 2008

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

Meditech Flow Selectors are Class IIa medical devices. CE marked in accordance with the European Medical Device regulations EU 2017 / 745. All items in this range are sold with a return to manufacturer warranty against manufacturing defects.

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website. Entonox is a registered Trademark of the Linde Group.

medITECH
MEDICAL GAS EQUIPMENT

Distributor: