

SPECIFICATION

Intended Gas:	Medical Oxygen (O ₂ , MA4, O ₂ /N ₂ O)
Inlet pressure:	4 bar (400 kPa)
Set Flow Rates in lpm:	Three flow rates options available: 1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm 10, 20, 30, 35, 40, 50, 55, 60, 70 lpm
Inlet Connections:	Schrader (BS 5682:2015) and other international connections
Outlet (Patient) Connections.	Firefree style barbed outlet, 9/16 threaded outlet, 3/8" BSP threaded outlet & 1/4" barb
Intended Life:	10 Years
Service Interval:	5 Years
Warranty Period:	5 Years
Complies With:	ISO 15002: 2008

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

Meditech Flow Selectors are Class IIa medical devices, CE marked in accordance with the European Medical Device Regulations EU 2017 / 745. All items in this range are sold with a return to manufacturer warranty against manufacturing defects.

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website. Entonox is a registered Trademark of the Linde Group.

medITECH
MEDICAL GAS EQUIPMENT

Distributor:

CE 2797

Doc 779M Issue A 27/04/21