

# medITECH

MEDICAL GAS EQUIPMENT

## Gas Specific Adaptor

Straight Bar Twin Adaptor - BS Schrader - 3/8" Nut + Liner



- Available with International Fittings.
- Lightweight construction.
- Available with 50mm or 100mm back bar.
- Service Free.
- 5 Year Warranty.

Excellent performance and robustness, coupled with a very lightweight design.

The Meditech Straight Bar Twin Adaptor is designed for use in all Hospital and Emergency Environments. This device has a Meditech 5 year warranty and is service free for its 10 year recommended life period.

Each Gas Specific Adaptor is available with the Outlet Types shown below. If the Outlet type you require is not shown, please contact us to discuss your requirements.

### OUTLET TYPES:



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS



CZECH



UNIFOR



### AVAILABLE GAS TYPES:

- O<sub>2</sub>
- MA4
- MA7
- O<sub>2</sub>/N<sub>2</sub>O
- N<sub>2</sub>O
- VACUUM

# SPECIFICATION

Intended Gas: .....	Medical Oxygen (O2, MA4, O2/N2O)
Inlet pressure: .....	4 bar (400 kPa)
Inlet Connections: .....	3/8" Nut + Nipple
Outlet Connections (4-BAR): .....	BS Schrader & Other interational Outlets.
Intended Life: .....	10 Years
Service Interval: .....	5 Years
Warranty Period: .....	5 Years

**TO PURCHASE OR FOR FURTHER DETAILS:**

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**  
**9 Fifth Avenue, Bluebridge Industrial Estate,**  
**Halstead, Essex, CO9 2SZ, England, UK.**  
**TEL: +44 (0) 1787 479 475**  
**FAX: +44 (0) 1787 477747**  
**EMAIL: sales@meditech.uk.com**  
**WEBSITE: www.meditech.uk.com**

All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Entonox is a registered Trademark of the Linde Group.

**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor: