

Gas Specific Adaptor

Twin Adaptor - BS Schrader - BS Probe



- **Available with International Fittings.**
- **Lightweight construction.**
- **Available with 50mm or 100mm back bar.**
- **Service Free.**
- **5 Year Warranty.**

Excellent performance and robustness, coupled with a very lightweight design.

The Meditech Straight Bar Twin Adaptor is designed for use in all Hospital and Emergency Environments. This device has a Meditech 5 year warranty and is service free for its 10 year recommended life period.

Each Gas Specific Adaptor is available with the Probe Types and / or the Outlet Types shown below. If either the Probe or Outlet type you require is not shown, please contact us to discuss your requirements.

PROBE TYPES:



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS



CZECH



UNIFOR



OUTLET TYPES:



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS




CZECH



UNIFOR



THERAPY TYPES:



FIR TREE



9/16" THREADED
THERAPY
CONNECTOR



3/8" BSP
THREADED THERAPY
CONNECTOR



1/4" BARB

AVAILABLE GAS TYPES:

- O₂
- MA4
- MA7
- O₂/N₂O
- N₂O
- VACUUM

SPECIFICATION

Intended Gas:	Medical Oxygen (O ₂ , MA4, O ₂ /N ₂ O)
Inlet pressure:	4 bar (400 kPa)
Inlet Connections:	Schrader Prove (BS 5682:2015) and other international connections
Outlet Connections.	BS Schrader Check Unit and other international outlets.
Intended Life:	10 Years
Service Interval:	5 Years
Warranty Period:	5 Years

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.**

TEL: +44 (0) 1787 479 475

FAX: +44 (0) 1787 477747

EMAIL: sales@meditech.uk.com

WEBSITE: www.meditech.uk.com

All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Entonox is a registered Trademark of the Linde Group.

medITECH
MEDICAL GAS EQUIPMENT

Distributor: