

## Pin Index Fixed Yolk Filling Adaptor (1/4" BSP)

- Adaptor for filling Pin Index Cylinders.
- Variety of thread types available - please specify.
- Available for any Medical Gas type.



### SPECIFICATION

Thread Type: ..... 1/4" BSP (other thread types available upon request)

Gas Type: ..... Oxygen (other gas types available upon request)

# medITECH

MEDICAL GAS EQUIPMENT

**TO PURCHASE OR FOR FURTHER DETAILS:**

Contact us by email or telephone for more details,  
including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**  
**9 Fifth Avenue, Bluebridge Industrial Estate,**  
**Halstead, Essex, CO9 2SZ, England, UK.**  
**TEL: +44 (0) 1787 479 475**  
**FAX: +44 (0) 1787 477747**  
**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**  
**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website. Entonox is a registered Trademark of the Linde Group.

Distributor: