

medITECH

MEDICAL GAS EQUIPMENT

Nitrous Oxide (N₂O) Regulator

Pin Index, 1 X BS Schrader, NON Therapy



699-0014-10 (Brass)

- Available with any international high pressure outlet fitting
- Luminous gauge face for night use
- 5-Year warranty
- 5-Year service interval

Meditech offers a wide range of high-pressure regulators for medical use and within this is our Nitrous Oxide models. They are available for use with any Outlet / Inlet Fitting needed to meet local requirements / specifications.

SPECIFICATION

Rated Inlet Pressure:	200 bar
Rated Outlet Pressure:	4 bar (400 kPa)
Storage temperature range:.....	-40°C to +60°C
Operating temperature range:	-20°C to +60°C
Intended Life:	10 Years
Service Interval:	5 Years
Warranty Period:	5 Years
Inlet Type:	Pin Index *Others available upon request*
Outlet Type:	BS Schrader *Others available upon request*

medITECH

MEDICAL GAS EQUIPMENT

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex CO9 2SZ, England, UK.**

Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747

E-mail: sales@meditech.uk.com • www.meditech.uk.com

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

Meditech Medical devices are CE marked in accordance with the Medical Devices Directive in the EU. All items in this range are sold with a return to manufacturer warranty against manufacturing defects.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd. terms and conditions of sale apply, available on application or from our website.

Distributor: