


## 2L Medical Oxygen Cylinder

BNOS Meditech has a Wholesalers Dealers License, allowing the sale and distribution of Medical Oxygen. Authorisation Holder Number: WDA (H) 44056



### SPECIFICATION

Filling Pressure:.....	200 BAR
Test Pressure:.....	300 BAR
Material:.....	Aluminum Alloy 6061
Diameter:.....	~114mm
Weight(filled with valve):.....	~4Kg
Weight(Empty with valve):.....	~2.8Kg
Height:.....	~ 444.5mm
Colour Code:.....	White* - EU Standard
Standard:.....	BS EN ISO 7866:2012
Valve: .....	Pin Index or Bullnose**
Volume (L):.....	400 Litres
Thread: .....	17E Tapered

\* Other colours available if required

\*\* Other valve types available if required

# medITECH

MEDICAL GAS EQUIPMENT

**TO PURCHASE OR FOR FURTHER DETAILS:**

Contact us by email or telephone for more details,  
including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**

**9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex, CO9 2SZ, England, UK.**

**TEL: +44 (0) 1787 479 475**

**FAX: +44 (0) 1787 477747**

**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**

**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.