

40L Medical Oxygen Cylinder

BNOS Meditech has a Wholesalers Dealers License allowing the sale and distribution of Medical Oxygen.
Authorisation Holder Number: WDA (H) 44056


Valve Types:



Pin Index



Bullose



SPECIFICATION

Filling Pressure:.....	200 BAR
Test Pressure:.....	300 BAR
Material:.....	Steel
Diameter:.....	~ 241mm
Weight(filled with valve):.....	~40Kg
Weight(Empty with valve):	~36.2Kg
Height:.....	~ 1295mm
Colour Code:.....	White* - EU Standard
Standard:.....	BS EN ISO 7866:2012
Valve:	Pin Index or Bullnose**
Volume (L):.....	8000 Litres
Thread:	25E Tapered

* Other colours available if required

** Other valve types available if required

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.