


## 4 Way Flow Selector On Stand Mount

BS Schrader Probe on 2M Supply Hose



- **No Calibration required EVER**
- **Lightweight construction**
- **Available with International Fittings**
- **7 Year Warranty**
- **7 Year service Interval**
- **Mass casualty device**


Excellent performance and robustness, coupled with a very lightweight design. Manufactured from a strong aluminium stand mount. Supplied with 2M supply hose to connect 4 BAR Oxygen source (longer hose length available upon Request).

For use in all Hospital and Emergency Environments, the 4 Way Flow Selector is designed to be plugged into Medical Oxygen Wall Outlets.


The Meditech 4 Way Flow Selector has been designed to deliver precision flow rates, without the need for any calibration. It has a standard 2M hose length (other hose lengths available upon request.)

Supplied with Push Fit "fir-tree" for direct attachment of Oxygen Therapy accessories.

### PROBE TYPES:



### THERAPY TYPES:



#### Available gas Types:

- O<sub>2</sub>
- MA4
- ENTONOX

# SPECIFICATION

Intended Gas: .....	Medical Oxygen (O <sub>2</sub> )
Inlet pressure: .....	4 bar (400 kPa)
Set Flow Rates in lpm: .....	Three flow rates options available: 1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm 10, 20, 30, 35, 40, 50, 55, 60, 70 lpm
Inlet Connections: .....	Schrader (BS 5682:2015) and other international connections
Outlet (Patient) Connections. ....	Firtree style barbed outlet.
Intended Life: .....	10 Years
Service Interval: .....	7 Years
Warranty Period: .....	7 Years
Complies With: .....	ISO 15002: 2008
Weight including hose:.....	1.6Kg
Dimensions: .....	Height: 15CM x 12CM (excluding hose)

**TO PURCHASE OR FOR FURTHER DETAILS:**  
**Contact us by email or telephone for more details,**  
**including specific part numbers, prices, etc.**

B.N.O.S. Meditech Ltd.  
9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex CO9 2SZ, England, UK.  
Tel: +44 (0) 1787 479475 • Fax: +44 (0) 1787 477747  
E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com) • [www.meditech.uk.com](http://www.meditech.uk.com)

Meditech Flow Selectors are CE marked in accordance with the Medical Devices Directive in the EU and UKCA to the Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website. Entonox® is a trade mark of the Linde Group

**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor: