



Non-Magnetic Medical Gas Regulator

- **7 Year Warranty**
- **Available with international fittings**
- **For use in environments that require non-magnetic devices.**
- **Service free over product lifetime**
- **Light-weight construction**
- **0-15 & 0-25lpm available**

Meditech offers a wide range of Regulators for medical use and within this is our Non magnetic models. Used in environments such as MRI Scans.

The Photo-Luminescent Gauge allows for easy reading, even in dimly lit environments.

Low cost of ownership and maintenance free design.



SPECIFICATION

Inlet Pressure:	8-200 bar (800-20000 kPa)
Set Flow Rates in lpm:.....	0-15 or 0-25lpm
Inlet Connections:	Pin Index, Bullnose & CGA-540
Outlet (patient) connection:	Fir tree barb or 9/16" threaded fitting (for use with humidifiers)
Intended Life:	10 Years
Service Interval:.....	Service free over product lifetime
Warranty Period:	7 Years
Complies with:	ISO 10524-1-2019

medITECH

MEDICAL GAS EQUIPMENT

B.N.O.S. Meditech Ltd.

9 Fifth Avenue, Bluebridge Industrial Estate,

Halstead, Essex CO9 2SZ, England, UK.

Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747

E-mail: sales@meditech.uk.com • www.meditech.uk.com

TO PURCHASE OR FOR FURTHER DETAILS:

**Contact us by email or telephone for more details,
including specific part numbers, prices, etc.**

B.N.O.S. Meditech Ltd.

9 Fifth Avenue, Bluebridge Industrial Estate,

Halstead, Essex CO9 2SZ, England, UK.

Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747

E-mail: sales@meditech.uk.com • www.meditech.uk.com

Meditech Regulators are CE marked in accordance with the Medical Devices Directive in the EU and UKCA to the Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor: