

CTS Splitter Unit

Combined Therapy Suction

The Meditech CTS Splitter Unit is a fully integrated oxygen and suction supply device which is powered from an gas cylinder. It is predominantly used on patient transfer or resuscitation trolleys in a hospital setting or for use where piped oxygen or vacuum are not available.

It can be supplied with numerous mounting options and hose lengths and is compatible with all types of collection jars / containers.



Available with Click Stop Flow Selector.

Available with Bobbin Style Flow Selector.

SPECIFICATION

Intended Gas:	Medical Oxygen (O ₂) or Medical Air (MA4)
Inlet pressure:	4 bar (400 kPa)
Max suction range:.....	40 lpm
Oxygen Flow Rates:	0-15 lpm
Inlet Connections:	Schrader (BS 5682:2015) and other international connections
Intended Life:	10 Years
Service Interval:	7 years

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

Distributor:

All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Doc 1014M Issue A 16/04/2025