

# meditech

## MEDICAL GAS EQUIPMENT

### Top Entry Regulators



**Bullnose Inlet**



**Pin Index Inlet**

- 7 Year Warranty
- Available with international Inlet & Outlet fittings
- Service free over product lifetime
- Light-weight construction
- Luminous Gauge face for night use
- 0-15,0-25 & 0-70lpm available

Meditech offers a wide range of high-pressure regulators for medical use and within this is our Top Entry Models. They are available for use with any outlet / inlet fitting needed to meet local requirements / specifications.

## SPECIFICATION

Inlet Pressure: .....	200 bar
Outlet Pressure: .....	4-bar (400 Kpa)
Inlet Connections: .....	Pin Index or Bullnose Others available upon request
Outlet (patient) connection: .....	Fir Tree or 9/16" threaded fitting Others available upon request
Intended Life: .....	10 Years
Service Interval:.....	Service free over product lifetime
Warranty Period: .....	7 Years

# medITECH

MEDICAL GAS EQUIPMENT

**B.N.O.S. Meditech Ltd.**

**9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex CO9 2SZ, England, UK.**

**Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747**

**E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com) • [www.meditech.uk.com](http://www.meditech.uk.com)**



Meditech Regulators are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. Meditech is the exclusive UK distributor for Flow Meter Italy and their Bobbin Style Flow Meter Range.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor: