



Pin Index Oxygen Regulator Side Entry Incl. Bobbin Flow Meter

- 7 Year Warranty
- Available with international fittings
- Service free over product lifetime
- Light-weight construction
- 0-1, 0-4, 0-15,0-25 & 0-70lpm available

Meditech offers a wide range of high-pressure regulators for medical use and within this are our Side Entry models. They are available for use with any outlet / inlet fitting needed to meet local requirements / specifications.

SPECIFICATION

| | |
|------------------------------------|---|
| Inlet Pressure: | 200 bar |
| Set Flow Rates in lpm:..... | 0-1lpm, 0-4 lpm, 0-15,0-25 & 0-70lpm available |
| Inlet Connections: | Pin Index Others available upon request |
| Outlet (patient) connection: | Fir Tree or 9/16" threaded fitting Others available upon request |
| Intended Life: | 10 Years |
| Service Interval:..... | Service free over product lifetime |
| Warranty Period: | 7 Years |



B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex CO9 2SZ, England, UK.**

Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747

E-mail: sales@meditech.uk.com • www.meditech.uk.com



Meditech Regulators are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. Meditech is the exclusive UK distributor for Flow Meter Italy and their Bobbin Style Flow Meter Range.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor: