

Bobbin Style Flow Selector

With Splitter - Direct Probe, 1 X 4-BAR Outlet



- **No Calibration required EVER**
- **Lightweight construction**
- **Available with all International Fittings**
- **2 Year Warranty**
- **10 year service free**
- **Available in Oxygen or Medical Air**
- **Supplied with a 9/16" or 3/8" BSP threaded outlet to screw on fir tree connector or humidifier.**

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Bobbin Style Flow Selector with Splitter is designed to be plugged into Medical Gas Wall Outlets.

The Meditech Bobbin Style Flow Selector with Splitter has been designed to deliver precision flow rates, without the need for any calibration.


Can be supplied with 5 flow rates:
0-1lpm, 0-4lpm, 0-15lpm, 0-30lpm & 0-50lpm

Our Bobbin Style Flow Selector with Splitter is available with the Probe types and Outlet types shown below. If either the Probe type or Outlet type you require is not shown, please contact us to discuss your requirements.

PROBE TYPE



OUTLET TYPES



SPECIFICATION

Intended Gas:	Medical Oxygen (O2), Medical Air (MA4)
Inlet pressure:	4 bar (400 kPa)
Set Flow Rates in lpm:	Three flow rates options available: 0-1lpm, 0-4lpm, 0-15lpm, 0-30lpm & 0-50lpm
Inlet Connections:	BS probe or any other international probe types.
Outlet (Bobbin) Connections.....	9/16" threaded outlet, 3/8" BSP threaded outlet, supplied with screw on therapy barb
Outlet (Splitter) Connections.....	BS Schrader or any other international probe types.
Intended Life:	10 Years
Service Interval:	10 year service free
Warranty Period:	2 Years
Complies With:	ISO 15002: 2008

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.**

TEL: +44 (0) 1787 479 475

FAX: +44 (0) 1787 477747

EMAIL: sales@meditech.uk.com

WEBSITE: www.meditech.uk.com



Meditech Gas Specific Adaptors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

medITECH
MEDICAL GAS EQUIPMENT

Distributor:

Doc 1027M Issue A 26/06/2025