mediTECH
MEDICAL GAS EQUIPMENT

2L Medical Oxygen Cylinder



B.N.O.S Meditech is fully authorised to sell and distribute Medical Oxygen under Wholesale Distribution Authorisation Number: WDA (H) 44056.

Valve Types:



Pin Index



Bullnose

SPECIFICATION

Filling Pressure:	200 BAR
Test Pressure:	300 BAR
Material:	Steel
Height:	~410mm
Diameter:	~113mm
Weight (empty with valve):	~5kg
Weight (filled with valve):	~6kg
Colour Code:	White* - EU Standard
Standard:	BS EN ISO 7866:2012
Valve:	Pin Index or Bullnose**
Volume (L):	400 Litres
Thread type:	17E Tapered or 25E Tapered

* Other colours available if required
** Other valve types available if required

mediTECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:
Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.